



Relaxatio™



Treatment - simple, fast and safe

1. Make sure the Relaxatio is fully charged before use (green light LED indicator).
2. After you turn on the device, it will start heating up slowly (no need to wait).
Place the glass probe on your skin and the red light will automatically activate.
3. Apply light pressure while slowly moving the Relaxatio in small circular motions.
4. Treat the painful areas: lower back, neck, shoulder, etc.
5. It is recommended to use the Relaxatio with a small amount of pain relief cream or gel (Red Algae Infused X Pain Relief Gel) on the treated areas during each session and to massage it gently in circular motions. Any cream you use with the treatment will make it more effective.
6. Clean the Relaxatio after every use with a cloth, cotton ball or paper towel.